TECH. TINCTUREOF GREEN SOAP

TECHNICAL DATA SHEET

GREEN SOAP IN 30% ETHANOL SOLUTION + LAVENDER OIL

DESCRIPTION:

Per USP & National Formulary guidelines, Sanitek uses only natural oils derived from vegetables and is free from coconut or coconut palm kernel oils. Technical Tincture of Green Soap contains no animal products, dyes, or fillers.

Technical Tincture of Green Soap is a mild, effective solution of green soap (aka soft soap, sapo mollis, sapo viridis) in ethanol for general cleaning of medical instruments and utensils and skin washing.

Natural and biodegradable, Technical Tincture of Green Soap contains glycerin that helps to restore moisture to skin and minimize the drying effects of alcohol.

INGREDIENTS:

Green Soap, 30% SD Ethanol, Glycerin, Lavender Oil.

FOR USE IN:

Traditional bulk-fill soap dispensers. Always dilute 50% with soft or de-ionized water when using in foaming dispensers.

USAGE:

For External Use Only.

Dilute up to 9:1 for general hand & skin washing and sanitation

May also be used for cleaning instruments and utensils.

HEALTH & SAFETY:

WARNING: FLAMMABLE LIQUID & VAPOR

DO NOT TAKE INTERNALLY • HARMFUL IF SWAL-LOWED OR INHALED

CAUSES IRRITATION TO EYES AND RESPIRATORY TRACT.

HEALTH & SAFETY: Cont.

AFFECTS CENTRAL NERVOUS SYSTEM.

MAY BE HARMFUL IF ABSORBED THROUGH SKIN

MAY CAUSE IRRITATION TO SKIN.

See SDS for More Health & Safety Information

For Industrial and Institutional Use Only

KEEP OUT OF THE REACH OF CHILDREN

Do Not Store in Open or Unmarked Containers

PACKAGING:

12/16oz Case, 4/1-Gallon Case, 55-Gallon Drum, 260-Gallon Tote

ITEM NO.: 420

NATIONAL FORMULARY:

Green Soap —

The Vegetable Oil	380.0 g
Oleic Acid	20.0 g
KOH, 85% active	91.7 g
Glycerin	<u>50.0 ml</u>
Purified Water to make	1,000 g

Green Soap Tincture —

Green Soap	650.0 g
Suitable essential oil(s) Alcohol	316.0 ml
Purified water to make	1,000 ml

Sanitek manufactures these soaps per above recipes in closed, pressure vessels systems operated under heat & pressure for unsurpassed quality.

Made in USA Rev. 03/14-1



3959 Goodwin Ave. Los Angeles, CA 90039

www.sanitek.com Phone: (323) 245-6781 Fax: (818) 242-1071 email: info@sanitek.com

This document has been prepared using data from sources considered technically reliable. It does not constitute a warranty, express or implied, as to the accuracy of the information contained herein. User is to determine suitability of this product as actual use conditions are beyond control of the seller. User is responsible to evaluate all available information when using product and to comply with all Federal, State and Local laws and regulations.